

ISSUE 35 - Published June 2025

£4.95

Animal Therapy

MAGAZINE

*Elevating education in animal
health & welfare*

BETA International

Sponsored by Neue Schule, returns reimagined

OMEGA-3S & JOINT HEALTH
Why Green-Lipped Mussel
Deserves a Closer Look

Unleashing
**HEALING
POWER:**

The **DEEP OSCILLATION®**
Revolution for
Horses & Hounds

Equicantis Business Live:
*A Day to Think Bigger,
Bolder, and Braver in
Your Therapy Business*

Tiny Tim

A HEARTWARMING
TALE OF RESILIENCE

Smart TECAR in **Science,**
Veterinary Medicine: *Not Magic*

Smart TECAR in Veterinary Medicine: *Science, Not Magic*

In the world of veterinary rehabilitation and physiotherapy, technology is evolving rapidly. Among the innovative tools now enhancing treatment outcomes is Smart TECAR therapy, also known as radiofrequency therapy, a modality that's revolutionising how we approach healing in animals. While some may mistakenly view this technology as "miraculous," the truth is far more grounded: Smart TECAR is not magic. It's science, and when used with knowledge and precision, it can dramatically optimise healing across a wide range of animal patients.

What is Smart TECAR (Radiofrequency) Therapy?

TECAR stands for "Transfer of Energy Capacitive and Resistive." It is a form of radiofrequency therapy that uses high-frequency electromagnetic energy to stimulate tissue repair, reduce pain, and promote circulation. The "Smart" component refers to enhanced features like real-time monitoring, adjustable power output, and feedback-driven treatment, tools that allow practitioners to tailor therapy to each patient with remarkable accuracy.


Smart TECAR therapy doesn't just warm the surface; it penetrates deep into tissues depending on the chosen mode (capacitive for superficial structures and resistive for deeper ones like joints and tendons). The result is accelerated metabolism, reduced inflammation, and faster healing, supported by solid physiological principles, not fantasy or placebo. It's a seamless blend of physics, bioelectricity, and biology working in sync.

Power in Watts: The Unsung Hero of Treatment Feedback

A key factor that separates Smart TECAR (radiofrequency) therapy from other tools is the ability to monitor power output in real time, measured in watts. This might seem like a technical detail, but it's a game changer in clinical practice.

Each animal, whether it's a racehorse, a German Shepherd, or a guinea pig, has different tissue

composition, hydration levels, and therapeutic needs. Applying a generic setting ignores these variables and risks inefficiency or discomfort. But by observing wattage during treatment, clinicians can fine-tune energy delivery on the spot, ensuring the optimal therapeutic dose is being delivered.



This is not guesswork. This is data-driven treatment in action. Rather than simply relying on time or subjective feedback, practitioners now have a quantifiable metric to guide energy application, making every session more targeted and effective.

Versatility from Stable to Surgery Table

One of the most impressive aspects of Smart TECAR radiofrequency therapy is its adaptability across species. In equine therapy, it supports muscle recovery, joint healing, and even injury prevention. For dogs and cats, it's a valuable tool for managing osteoarthritis, post-surgical rehab, disc disease, and more.

Even small animals and exotics, often overlooked in rehab modalities, can

benefit, as long as the clinician carefully adjusts parameters. That's where wattage feedback becomes crucial. It enables safe, species-specific therapy from the smallest pets to the largest livestock.

This universality means one device can support a wide range of clinical scenarios, a smart investment for any veterinary physiotherapy practice.




About That Gel — No Magic Formula Required

One persistent myth in TECAR therapy is the supposed need for specialised or proprietary gels. The truth? Any standard water-based gel works just fine. The purpose of the gel is to facilitate effective contact and conductivity between the applicator and the patient's skin, not to add any pharmacological effect.

Smart TECAR therapy relies on the conductivity of water molecules to transmit radiofrequency energy into the body. If a gel is water-based and maintains smooth contact, it does its job. There is no “secret sauce” here, and no need for expensive branded gels. Focus your investment where it counts, on technique, equipment, and training.

A Fusion of Clinical Insight and Scientific Precision

Smart TECAR, or radiofrequency therapy, represents the next evolution of evidence-based, non-invasive treatment in veterinary medicine. But it is not a “magic wand.” Its effectiveness lies in the



science of controlled energy transfer, precise parameter management, and biological response, guided by the clinician's knowledge of anatomy and pathology.

Monitoring wattage in real time gives a window into how the tissue is responding, letting the clinician adjust the energy output to suit each patient's size, species, and condition. That kind of feedback isn't just useful, it's essential for optimising outcomes and ensuring safety.

And when it comes to treatment gel? Keep it simple. Choose a basic water-based gel, ensure smooth application, and let the technology, not marketing, do the work.

Conclusion: Real Science, Real Results

Smart TECAR radiofrequency therapy is a powerful tool in the hands of skilled and knowledgeable veterinary professionals. Its impact is measurable, its mechanisms are rooted in science, and its results are increasingly undeniable. But success doesn't come from mystery or magic formulas. It comes from understanding, adjusting, and applying — using evidence-based protocols, careful monitoring, and species-specific expertise.

When used correctly, Smart TECAR becomes more than a machine, it becomes an extension of clinical insight, empowering better care for animals of all shapes and sizes.



therapévo
the power to heal



Gamma Class IV Laser

Veterinary SMART TECAR v2



the #1 choice for all electronic therapy devices

www.therapevo.co.uk

class IV laser – shockwave – pemf – radio frequency – ultrasound – cryotherapy – class 3b laser

